


8.20 PARTS LIST AND EXPLODED VIEW

Model 2522C-01			
Item	Qty.	Part No.	Description
9	1	638122	Pressure Gauge
11	1	633210	Vacuum Gauge
14	1	638954-502	Pressure Regulator Ass'y
15	1	638755-502	Vacuum Regulator Ass'y

Model 2522B-01 & 2522Z-01 (115V 60Hz)			
Item	Qty.	Part No.	Description
9	1	638122	Pressure Gauge
11	1	633210	Vacuum Gauge
14	1	638947-502	Pressure Regulator Ass'y
15	1	638755-502	Vacuum Regulator Ass'y

Model 2522C-02 & 2522Z-02 (230V 50Hz)			
Item	Qty.	Part No.	Description
9	1	638787	Pressure Gauge
11	1	638792	Vacuum Gauge
14	1	638947-502	Pressure Regulator Ass'y
15	1	638755-502	Vacuum Regulator Ass'y

Model 2534B-01 (230V 50Hz)			
Item	Qty.	Part No.	Description
9	1	638122	Pressure Gauge
11	1	633210	Vacuum Gauge
14	1	638954-502	Pressure Regulator Ass'y
15	1	638755-502	Vacuum Regulator Ass'y



Item	Qty.	Part No.	Description	Service Kit 2522K-03	Seal Kit 2522K-04	Foot Kit 2500K-04
1	1	SK61722	Connecting Rod Ass'y	1	-	-
2	1	Note1	Screw, Handle	2	-	-
3	1	Note1	Screw, Head	4	-	-
4	1	Note1	O-ring, Head	1	1	-
5	1	621276-504	Valve Plate Ass'y	1	-	-
6	1	Note1	O-ring, Cylinder	1	1	-
7A	1	Note1	Jar - Plastic	1	1	-
7B	1	Note1	Jar - Plastic	1	1	-
8	1	Note1	Filter Element	1	1	-
9	1	Note1	Ball	1	1	-
12	1	Note1	Gaslet	2	2	-
13	1	Note1	Rubber Suction Feet	-	-	4
16	1	641020	Replacement Filters	3	3	-
17	1	Note1	O-ring, Cylinder	1	1	-

Model 2522C-10 (230V 60Hz)			
Item	Qty.	Part No.	Description
5	1	621818-502	Valve Plate Ass'y
9	1	638122	Pressure Gauge
11	1	638710	Vacuum Gauge
14	1	638947-502	Pressure Reg. Ass'y
15	1	638755-502	Vacuum Reg. Ass'y

Note 1: These parts are not available separately, but are supplied in service kits.