



1422A Vacuum Regulator




- **Chemical Resistant**
- **No mercury**
- **Wetted Parts are 100% Stainless Steel.**
- **Compact and Portable.**
- **Gauge does not require electricity or batteries**
- **Gauge is very economical**

The Welch Vacuum Regulator Gauge gives you complete control over your vacuum process. Vacuum control is easily done by adjusting the Vacuum Regulator knob to the desired pressure setting on the vacuum Gauge.

The Welch Vacuum Regulator is a small, portable unit, which is ideal for many applications requiring Vacuum Regulation.

When your vacuum source, (pump, house vacuum), is turned on you will have the complete capability of controlling your vacuum between 0 – 30”hg on the Vacuum Gauge.

Typical vacuum Regulator Setup



Vacuum Regulator Gauge is designed to accurately maintain the set vacuum pressure desired. All the wetted parts are Stainless Steel to insure chemical resistance and longevity.

Specifications and Ordering Information

Model		1422A	
Analog Gauge in.Hg.	0 - 30		
Wetted Parts	Stainless Steel		
Vacuum Regulator Valve	Adjustable		
Connection Hose Barb in.			
Inlet in.(mm)	1/4		
Outlet in. (mm)	1/4		
Overall Diementions			
L in	5.25		
W in	4.25		
H in	3.5		
Ship Weight, lbs (kg)	4 (1.8)		



7301 N. Central Avenue, Skokie, Illinois 60077
 Phone: 847-676-8800
 Fax: (Technical) 847-677-8606 (Ordering) 920-451-4397
 E-mail: welchvacuum@thomasind.com
www.welchvacuum.com