


### 8.14 Parts List & Exploded View for Dry Vacuum System





ITEM	QTY	PART NUMBER	DESCRIPTION	SERVICE KIT 2025K-06*
1	4	See Note 1	Hose Clamp 5/16"	4
2	2	See Note 1	Hose Clamp 1/8"	2
3	0.35 Ft	See Note 1	Tubing, Fluoroelastomer	0.35 Ft
4	1.60 Ft	See Note 1	Tubing, Blended	1.60 Ft
5	1.30 Ft	See Note 1	Tubing, PTFE Duality	1.30 Ft
6	1	See Note 1	Gasket	1
7	1	See Note 1	U-Clamp	1
8	1	726123	Vacuum Gauge 2" Dial (mbar / kPa)	-
8	1	726125	Vacuum Gauge 2" Dial (in Hg / cm hg)	-
9	2	1423B	Glass Jarr Assembly	-
10	2	1415B	Replacement Glass Jarr	-

Note: 1. These parts are not available seperately,  
but they are supplied in the 2025K-06 Service Kit.


\* for Models 2025, 2026, 2027 & 2028 with S/N > ECO91300

8.15 Parts List & Exploded View for Repair Kits



Repair Kit 2037K-01  
(for Models: 2025)

This kit provides enough parts to service one side/head of this vacuum pump.



Repair Kit 2047K-01  
(for Model: 2026, 2027 & 2028)

This kit provides enough parts to service one side/head of this vacuum pump.