



PARTS LIST FOR MODEL 1405

ITEM	QTY	P/N	A*	B*	DESCRIPTION	ITEM	QTY	P/N	A*	B*	DESCRIPTION
1	3	2-00-2706			SCREW PAN HD 10-32 x 3/8" LG.	31	1	41-1736			GAS BALLAST VALVE
2	3	2-00-6712			SCREW SOCKET HD 10-32 x 3/4" LG.	32	1	41-1742	1	1	O-RING #011 BUNA-N
3	10	2-01-0114			SCREW HEX HD 1/4-20 x 7/8" LG.	33	3	41-1766	3	3	THRUST DISK
4	6	2-01-0316			SCREW HEX HD 5/16-18 x 1" LG.	34	1	41-2045			SMALL END PLATE, FINISHED
5	6	2-01-0356			SCREW HEX HD 5/16-18 x 3-1/2" LG.	35	1	41-2074			PULLEY w/ SET SCREW
6	1	2-01-5106			SCREW SOCKET HD 1/4-20 x 3/8" LG.	36	1	41-2092			ROTOR, INTAKE
7	4	2-01-6112			SCREW SOCKET HD 1/4-20 x 3/4" LG.	37	1	41-2093			SHAFT
8	2	2-01-6320			SCREW SOCKET HD 5/16-18 x 1-1/4 LG.	38	1	41-2153	1	1	EXHAUST VALVE, UPPER
9	1	2-01-8128			SETSCREW 1/4-20 x 1-3/4" LG.	39	1	41-2154	1	1	EXHAUST VALVE, LOWER
10	3	2-31-0112			HEX NUT 1/4-20, STEEL	40	2	41-2158	2	2	VALVE SPRING
11	1	2-36-9900			HEX LOCK NUT 1/4-20	41	1	41-2159	1	1	VALVE BLOCK SPACER
12	3	2-61-0000			WASHER, STEEL 3/16	42	1	41-2160			BAFFLE ASSEMBLY
13	10	2-63-0193			WASHER, LOCK 1/4" MEDIUM	43	1	41-2175			SETSCREW 1/4-20 x 2.25" LG
14	1	4-06-0754	1	1	ROTOR CLIP	44	6	41-2363			WASHER 5/16 ID x 1/2 OD
15	1	4-40-1200	1	1	STEEL BALL 3/8" DIA.	45	1	41-2757			LOCK NUT
16	1	41-0234	1	1	GASKET, INTAKE CHAMBER COVER	46	1	61-8101			OIL CASE ASSEMBLY
17	1	41-0403	1	1	GASKET, OIL CASE	47	1	61-8107A			SHAFT SEAL ASSEMBLY
18	1	41-0508	1	1	WASHER, THRUST	48	1	61-8108D			EXHAUST RING
19	1	41-0612			DUST CAP	49	1	61-8250C			CENTER PLATE
20	2	41-0613	2	2	WOODRUFF KEY #404	50	1	61-8309B			PRESSURE RELEASE TUBE #1
21	2	41-0624	2	2	WOODRUFF KEY	51	1	61-8310A			PRESSURE RELEASE TUBE #2
22	1	41-0643			GASKET, SEAL	52	1	61-8313A			VALVE BLOCK
23	1	41-0660	1	1	INTAKE SCREEN	53	2	61-8314A	2	2	VALVE
24	1	41-0672			SHAFT END PLATE	54	1	61-8329A			BRACKET, INTAKE SCREEN
25	2	41-0959	2	2	VANE, LARGE	55	6	61-8456A	6	6	WASHER 5/16, COPPER
26	1	41-0992	1	1	PRESSURE RELEASE VALVE	56	1	61-8641			HOSE ADAPTER (7/16 HOSE I.D.)
27	6	41-1304	6	6	VANE SPRING	57	1	61-8651			INTAKE RING
28	2	41-1487	2		VANE, SMALL	58	1	61-8668			LARGE END PLATE
29	1	41-1518			ROTOR, EXHAUST	59	1	61-8671			COVER, FINISHED
30	3	41-1685	2	2	SPRING HOLDER	60	1	61-9284A	1	1	O-RING #214 VITON
						62	3	62-0259			SCREW #6-1/4, SELF-TAPPING
								1-99-5282	1	1	SEALANT, 6 cc
								1401D	1	1	MECHANICAL SEAL OPTION**

A* : MAJOR REPAIR KIT CAT# 1405K-06

B* : MINOR REPAIR KIT CAT# 1405K-05

** 1401D is mechanical seal includes mechanical seal, gasket and 3 screws