



Model 2012

67-1610  
Rev. 08 03/06

## PARTS LIST FOR MODEL 2012

ITEM	P/N	QTY	A*	B*	DESCRIPTION	ITEM	P/N	QTY	A*	B*	DESCRIPTION
1	2-00-6708	2			SCREW SOCKET HD 10-32 x 5/8"	20	61-3767	1			MOMENTUM WEIGHT
2	2-30-5708	2			NUT HEX 10-32 x 3/8 x 1/8	21	61-3769	1			HANDLE
3	2-62-0793	2			WASHER LOCK #10 MED.	22	61-3917	1			CONNECTING ROD ASSEMBLY
4	2-69-9901	4	1		LOCK WASHER 1/4 S.S.	23	62-0174	4			SCREW SOCKETHEAD M4 x 18
5	2-71-0040	6			WASHER BRASS FLAT #10L	24	62-0178	4			SCREW SOCKETHEAD M5x50
6	61-3059	1	1		DIAPHRAGM MOLD	25	62-0189	1			SET SCREW M5 x 8
7	61-3061	2	2		VALVE 8011	26	62-0195	4			SCREW SELF TAP 2.9x1x6
8	61-3093	1			HEAD PLATE	27	62-0200	24			BELLEVILLE WASHER 10x5
9	61-3102	2			DIAPHRAGM SHIM	28	62-0202	1			DOWEL PIN 3mm x 12mm
10	61-3479	1			DIAPHRAGM CHAMBER	29	62-1000	4			SOCKET HEAD M4x12
11	61-3480	1			INTAKE HEAD	30	62-1097	4			ACORN NUT #8-32
12	61-3482	1			DIAPHRAGM SUPPORT PLATE	31	62-1098	4			SCREW WASHERHEAD 8-32 x 1
13*	61-3518	1			MOTOR, 1/5 H.P. 230V 50Hz	32	62-1191	4	1		SCREW SOCKETHEAD M6 x 15
14	61-3553	1			BASE PLATE	33	62-1205	1			SET SCREW M10 x 10
15	61-3557	1			HOUSING LID, WHITE	34	66-0140	1			HOSE BARB 1/4 NPT x 3/8 HOSE
16	61-3560	1			HOUSING, FINISHED	35	66-0151	1			SILENCER, EXHAUST 1/4" NPT
17*	61-3562	1			MOTOR, 1/5 H.P. 115V 60Hz	36	66-0160	4			SUCTION CUP
18*	61-3563	1			MOTOR, 1/5 H.P. 100V 50/60Hz	37	66-0166	2	2		O-RING #021, FLUOROELASTOMER
19	61-3586	1			COVER PLATE	38	66-0170	4			RUBBER STOPPER

A\*: MAJOR REPAIR KIT#2025K-01 (S/N>EG050000)

A\*: MAJOR REPAIR KIT#2022K-01 (S/N>EG040000)

B\*: MINOR REPAIR KIT# N/A

\*: Depending on Model Number

\*\*: Shipped loose