



**PARTS LIST
MODEL 2017 VACUUM PUMP**

ITEM	QTY	PC.NO.	A* DESCRIPTION	ITEM	QTY	PC.NO.	A* DESCRIPTION
01*	1	61-3562	Motor 115V 60 Hz	23	1	61-5047	Connector Barbed
01*	1	61-3518	Motor 220V 50 Hz	24	1	62-0119	O-Ring Tetraseal
01*	1	61-3563	Motor 100V 50/60 Hz	25	4	62-0195	Screw, Pan Head
02	4	62-1098	Hex Washer Head Screw	26	2	61-3479	Diaphragm Chamber
03	4	2-71-0040	Flat Brass Washer	27	2	61-3059	1 Diaphragm Mold
04	1	61-3553	Base, Machining	28	4	2-69-9901	Lock washer
05	4	66-0160	Suction Cup	31	2	62-0202	Pin Dowel
06	8	62-0178	Screw Socket Head	36	1	61-3258	Handle
07	48	62-0200	Belleville Washer	37	1	61-3511	Connecting Rod Assembly
08	2	61-3093	Head Plate	44	1	62-0190	Screw, Socket Head
09	4	62-1097	#8-32 Acorn Nut	45	1	61-3767	Momentum Weight
10	1	62-0189	Screw Set Socket Head	49	1	61-3840	Nut Assembly
11	4	61-3061	2 8011 Valve	50	2	62-0216	Screw Hex Head Cap
12	4	66-0166	2 O-Ring	51	2	62-1054	Locker Washer
13	1	66-0127	Tube 3/8" OD x 1/4" ID	52	2	62-1053	Screw Socket Head
14	1	61-3557	Lid Housing	53	1	61-3561	Handle Bracket
15	1	66-0131	Elbow Male Union	55	4	62-1191	Screw Socket Head
16	8	62-0174	Screw, Socket Head	56	2	41-2363	Washer
17	2	61-3482	Base Plate	57**	1	66-0140	Adapter
18	4	61-3102	Shim Diaphragm	58**	1	66-0126	Elbow 1/4NPTx1/4Male/Female
19	1	61-3560	Housing Casting	59**	1	66-0151	1/4" Exhaust Silencer
20	2	61-3480	Head Intake	60***	1	1412B	Silencer Kit
22	1	66-0125	Angle Right "T" natural			1-99-5126	1 Sealant
						62-0179	1 Screw

A*: MAJOR REPAIR KIT CAT # 2025K-01

*: Depending on the model number

**: Shipped Loose

***: Contains items 58 and 59 and may be substituted for item 57