



PARTS LIST FOR MODEL 2022

ITEM	P/N	QTY	A*	B*	DESCRIPTION	ITEM	P/N	QTY	A*	B*	DESCRIPTION
1	2-69-9901	4			LOCK WASHER 1/4 S.S.	24	62-0174	8			SCREW SOCKETHEAD M4 x 18
2	2-71-0040	4			WASHER BRASS FLAT #10L	25	62-0178	8			SCREW SOCKETHEAD M5x50
3	41-2363	2			WASHER 5/16 ID x 1/2 OD	26	62-0189	1			SET SCREW M5 x 8
4	61-3059	2	1		DIAPHRAGM MOLD	27	62-0190	1			SET SCREW M6 x 8
5	61-3061	4	2		VALVE 8011	28	62-0195	4			SCREW SELF TAP 2.9x1x6
6	61-3093	2			HEAD PLATE	29	62-0200	48			BELLEVILLE WASHER 10x5
7	61-3102	4			DIAPHRAGM SHIM	30	62-0202	2			DOWEL PIN 3mm x 12mm
8	61-3258	1			HANDLE	31	62-0216	2			SCREW HEXHEAD M8x16
9	61-3479	2			DIAPHRAGM CHAMBER	32	62-1044	1			O-RING, VITON 3mm ID x 1.5mm
10	61-3480	1			INTAKE HEAD	33	62-1053	2			SCREW SOCKETHEAD #10 x 5/8
11	61-3481	1			EXHAUST HEAD	34	62-1054	2			LOCK WASHER #10
12	61-3482	2			DIAPHRAGM SUPPORT PLATE	35	62-1097	4			ACORN NUT #8-32
13	61-3483	1			GAS BALLAST SCREW	36	62-1098	4			SCREW WASHERHEAD 8-32 x 1
14*	61-3518	1			MOTOR, 1/5 H.P. 230V 50Hz	37	62-1191	4			SCREW SOCKETHEAD M6 x 15
15	61-3553	1			BASE PLATE	38**	66-0126	1			ELBOW 1/4"NPT MALE TO FEMALE
16	61-3557	1			HOUSING LID, WHITE	39	66-0127	1			TEFLON TUBING
17	61-3560	1			HOUSING, FINISHED	40	66-0131	2			ELBOW UNION 1/4NPT x 3/8TUBE
18	61-3561	1			HANDLE BRACKET	41**	66-0140	1			HOSE BARB 1/4 NPT x 3/8 HOSE
19*	61-3562	1			MOTOR, 1/5 H.P. 115V 60Hz	42	66-0141	1			ELBOW 1/4 NPT x 3/8 BARB
20*	61-3563	1			MOTOR, 1/5 H.P. 100V 50/60Hz	43**	66-0151	1			SILENCER, EXHAUST 1/4" NPT
21	61-3766	1			CONNECTING ROD ASSEMBLY	44	66-0160	4			SUCTION CUP
22	61-3767	1			MOMENTUM WEIGHT	45	66-0166	4	2		O-RING #021, FLUOROELASTOMER
23	61-3840	1			NUT ASSEMBLY						

A*: MAJOR REPAIR KIT#2025K-01 (S/N>EG050000)
A*: MAJOR REPAIR KIT#2022K-01 (S/N>EG040000)
B*: MINOR REPAIR KIT# N/A
*: Depending on Model Number
**: Shipped loose