


6.14 Parts List and Exploded View for Pump Model 2561



Model 2561

Item Number	Qty.	Part No.	Description	Service Kit 2561K-03	Seal Kit 2561K-04	Feet Kit 2500K-04
1	2	Note 1	Piston Cup	2	-	-
2	2	Note 1	Cup Retainer Screw	2	-	-
3	2	Note 1	Cylinder Sleeve	2	-	-
4	2	Note 1	O-Ring, Gasket	2	2	-
5	6	Note 1	Head Screw	6	-	-
6	1	Note 1	Valve, Intake / Intake Stage	1	1	-
7	3	Note 1	Valve, Intake & Exhaust / Both Stages	3	3	-
8	2	Note 1	O-Ring, Cylinder Sleeve	2	2	-
10	4	Note 1	Rubber Suction Feet	-	-	4

Note 1: These parts are not available separately, but are supplied within the kits shown in last three columns.