



67-1358  
Rev. 2 11/02

Model 2567

**PARTS LIST  
MODEL 2567 VACUUM PUMP**

<b>ITEM</b>	<b>QTY</b>	<b>PART NO.</b>	<b>DESCRIPTION</b>	<b>SERVICE KIT 2567K-03</b>	<b>SERVICE SEAL KIT 2567K-04</b>	<b>SERVICE FEET KIT 2500K-04</b>
1	2	2567K-01	Connecting Rod Assembly	2	-	-
2	2	2567k-02	Valve Plate Assembly	2	-	-
3	2	62-1050	Cylinder O-Ring	2	2	-
4	2	62-1052	Valve Plate Gasket	2	2	-
5	1	1412C	Muffler	1	-	-
6	8	62-1051	Head Screw	8	-	-
7	4	66-0160	Feet, Suction Cup	-	-	4
8	2	62-1109	Retainer Hold Down Screw	-	2	-
9	2	62-1110	PTFE Piston Seal	-	2	-