


6.16 Parts List and Exploded View for Pump Model 2580



Model 2580

Item Number	Qty.	Part No.	Description	Service Kit 2585K-03	Seal Kit 2585K-04	Feet Kit 2500K-04
1	2	2585K-01	Connecting Rod Ass'y	2	-	-
2	2	Note 1	O-Ring Valve Plate	2	2	-
3	2	Note 1	O-Ring, Connector Tube	4	4	-
4	8	Note 1	Head Screw	8	-	-
5	4	Note 1	Rubber Suction Feet	-	-	4
6	2	2585K-02	Valve Plate Ass'y	2	-	-
7	2	Note 1	O-Ring Cylinder	2	2	-
8	1	1412C	Muffler	1	-	-

Note 1: These parts are not available separately, but are supplied within the kits shown in last three columns.