



Model 8111

67-1205
 Sheet 2 of 2
 Rev. 02 11/03

PARTS LIST FOR 8111

ITEM	QTY	PC. NO.	A*	DESCRIPTION
01	2	61-3403		Tube
02	1	61-3392		Tube Fitting
03	4	61-3393		Screw
04	2	61-3394		Valve Plate
05	2	61-3395		O-Ring
06	4	61-3396		Valve
07	8	61-3397		Screw
08	8	61-3398		Washer
09	2	61-3400		Pump Head
10	1	61-3124		Shaft
11	2	61-3401		Diaphragm 8111
12	2	61-3402		Shim, Diaphragm
13	1	66-0181		Fan
14	1	62-0198		Retaining Clip
15	2	-		Transition Plate

67-1205
 Sheet 1 of 2
 Rev. 2 11/03