

Aspiration System Model 2515

Replacement Jars for Aspiration Stations

Spare receiver kit (1475K-20) includes a 1.2 liter autoclavable jar, 6ft. of tubing and two hydrophobic filters.

Pump Model	CAT. No.
2511B-75, 2511C-75	1475K-20
2515B-75, 2515C-75	



2511

Repair Kits for Model 2515

ID	Description	Model 2515B-75 Kit CAT. No.	Models 2515C-75/2515C-76 Kit CAT. No.
1	Connecting Rod Assembly	----	----
2	Valve Plate Assembly	2515K-02	2515K-22

